



Product specifications

Name Amphetamine-HRP Antigen

Description Synthetic Amphetamine, HRP conjugated

Amino acid sequence -

Product host -

Product code 170035

Product formulation Liquid; Antigen-HRP conjugate

Product buffer solution 1:1 0.1 M Tris base / FBS, pH 7.0, 0.016% Gentamicin sulfate, 0.05% thimerosal

Reconstitution -

Shelf life and storage Shelf life unspecified, storage at 2–8 °C. Do not freeze.

Analyte description Amphetamine is a central nervous system stimulant drug used to treat ADHD and narcolepsy.

Product concentration -

Miscellaneous This product was previously available with a Catalog code P91-99-71H. It has been produced

using the same conjugation protocol.

Matched pair C03-99-72M-P (capture)