

Cannabinoid (THC)-HRP Antigen

Catalog number	170034
Description	Synthetic Cannabinoid, HRP conjugated
Amino acid sequence	-
Product host	-
Product formulation	Liquid; Antigen-HRP conjugate
Product buffer solution	1:1 0.1 M Tris base / FBS, pH 7.0, 0.0016% Gentamicin sulfate, 0.005% thimerosal
Reconstitution	-
Shelf life and storage	Shelf life unspecified, storage at 2–8 °C. Do not freeze.
Analyte description	Tetrahydrocannabinol (THC) is a psychoactive cannabinoid.
Product concentration	-
Matched pairs	100892, 100893, 100894, P01-99-21M-P, AG309 (capture)
Miscellaneous	This product was previously available with a Catalog code P91-99-21H. It has been produced using the same conjugation protocol.

**Legal disclaimer**