

## Product specifications

---

Name	Cannabinoid (delta-9-THC) Antibody
Specificity	Antibody primarily recognizes delta-9-THC with some reactivity to metabolites. Recommended for oral fluid testing applications.
Description	Protein A purified monoclonal mouse antibody, cultured <i>in vitro</i> under conditions free from animal-derived components
Product code	100893
Product buffer solution	0.05M Potassium Phosphate Buffer, pH 7.2 with 0.85% NaCl, 0.05% NaN <sub>3</sub> as a preservative
Shelf life and storage	12 months from manufacturing at 2–8 °C
Subclass	IgG <sub>1</sub>
Analyte description	Tetrahydrocannabinol (THC) is a psychoactive cannabinoid.

## Parameters tested on each lot

---

Product appearance	Liquid, may turn slightly opaque during storage
Product concentration	≥ 1.5 mg/ml
Immunoreactivity	Not Determined (N/D)
IEF Profile	N/D
Purity	≥ 90 %

## Kinetic parameters

---

Association rate constant	N/D
Dissociation rate constant	N/D
Affinity constant	N/D
Determination method	-
Determination antigen	-

### Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored and used according to the instructions. MedixMAB and MedixAntigens are trademarks of Medix Biochemica and may not be used or reproduced without Medix Biochemica's written permission.



Cross-reactivities	Does not react with Opiates, Benzoylcegonine, Amphetamines, Phencyclidine	
Epitope	N/D	
Pair recommendations	P81-99-21A (BSA conjugate as capture)	
	Please note that pair recommendations are based on results obtained by our laboratory. Equally good results may be obtained using other pairs and therefore these recommendations are only indicative.	
Platforms tested	ELISA, LF	
Antigens tested	Please see Pair recommendations above	
Product stability	TEMPERATURE, TIME	RESULT
	N/D	-
Miscellaneous	This product was previously available with Catalog codes P01-99-24M-P and AG310. It has been produced using the same original mouse monoclonal hybridoma cell line.	
	Immunogen: THC-BgG	
References	-	

**Legal disclaimer**

Medix Biochemica products meet their specifications if transported, stored and used according to the instructions. MedixMAB and MedixAntigens are trademarks of Medix Biochemica and may not be used or reproduced without Medix Biochemica's written permission.

