

Product specifications

Name	Tricyclic Antidepressant (TCA)-BSA Antigen
Description	Synthetic Tricyclic Antidepressant, BSA conjugated
Amino acid sequence	-
Product host	-
Product code	170030
Product formulation	Liquid; Antigen-BSA conjugate
Product buffer solution	Distilled water
Reconstitution	-
Shelf life and storage	Shelf life unspecified. Short term storage (≤ 2 weeks): 2–8 °C . Long term: -20 °C. Avoid repeated freezing and thawing.
Analyte description	Tricyclic antidepressants are a class of medications used primarily for the treatment of depression and other mood disorders.
Product concentration	Lot-specific
Miscellaneous	This product was previously available with Catalog codes T001A and CB500. It has been produced using the same conjugation protocol.
Matched pair	T002M or CB300 (detection).

Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored and used according to the instructions. MedixMAB and MedixAntigens are trademarks of Medix Biochemica and may not be used or reproduced without Medix Biochemica's written permission.

