

## Product specifications

---

Name	Fentanyl-BSA Antigen
Description	Synthetic Fentanyl, BSA conjugated
Amino acid sequence	-
Product host	-
Product code	170040
Product formulation	Liquid; Antigen-BSA conjugate
Product buffer solution	Distilled water
Reconstitution	-
Shelf life and storage	Shelf life unspecified. Short term storage ( $\leq 2$ weeks): 2–8 °C. Long term: -20 °C. Avoid repeated freezing and thawing.
Analyte description	Fentanyl is a synthetic opioid with 50–100 times the potency of morphine. It is used to treat severe pain. Fentanyl may be illicitly produced and mixed with other illegal drugs increasing the abuse and overdose risks.
Product concentration	Lot-specific
Miscellaneous	This product was previously available with a Catalog code P81-99-51A. It has been produced using the same conjugation protocol.
Matched pair	P01-99-53R-IF, P01-99-54R-IF, P01-99-55R-IF (detection)

**Legal disclaimer**

Medix Biochemica products meet their specifications if transported, stored and used according to the instructions. MedixMAB and MedixAntigens are trademarks of Medix Biochemica and may not be used or reproduced without Medix Biochemica's written permission.

