

Product specifications

Name	Zolpidem-BSA Antigen
Description	Synthetic Zolpidem, BSA conjugated
Amino acid sequence	-
Product host	-
Product code	170041
Product formulation	Liquid; Antigen-BSA conjugate
Product buffer solution	Distilled water
Reconstitution	-
Shelf life and storage	Shelf life unspecified. Short term storage (≤ 2 weeks): 2–8 °C. Long term: -20 °C. Avoid repeated freezing and thawing.
Analyte description	Zolpidem is a central nervous system depressant drug causing drowsiness, sedation and decreased psychomotor function. Zolpidem is used to treat insomnia and may result in next-day impairment.
Product concentration	Lot-specific
Miscellaneous	This product was previously available with a Catalog code P83-99-11A. It has been produced using the same conjugation protocol.
Matched pair	P03-99-11R-IF, P03-99-12R-IF, P03-99-13R-IF, JP046, JP047, JP048, HM311 or HM312 (detection)

Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored and used according to the instructions. MedixMAB and MedixAntigens are trademarks of Medix Biochemica and may not be used or reproduced without Medix Biochemica's written permission.

