

Product specifications

Name	Amphetamine-HRP Antigen
Description	Synthetic Amphetamine, HRP conjugated
Amino acid sequence	-
Product host	-
Product code	170035
Product formulation	Liquid; Antigen-HRP conjugate
Product buffer solution	Tris-HCl with enzyme stabilizers, pH 7.0
Reconstitution	-
Shelf life and storage	Shelf life unspecified, storage at 2–8 °C. Do not freeze.
Analyte description	Amphetamine is a central nervous system stimulant drug used to treat ADHD and narcolepsy.
Product concentration	-
Miscellaneous	This product was previously available with a Catalog code P91-99-71H. It has been produced using the same conjugation protocol.
Matched pair	C03-99-72M-P (capture)

Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored and used according to the instructions. MedixMAB and MedixAntigens are trademarks of Medix Biochemica and may not be used or reproduced without Medix Biochemica's written permission.

