

Product specifications

Name	Clonazepam-BSA Antigen
Description	Synthetic Clonazepam, BSA conjugated
Amino acid sequence	-
Product host	-
Product code	170064
Product formulation	Liquid; Antigen-BSA conjugate
Product buffer solution	Distilled water
Reconstitution	-
Shelf life and storage	Shelf life unspecified. Short term storage (≤ 2 weeks): 2–8 °C. Long term: -20 °C. Avoid repeated freezing and thawing.
Analyte description	Clonazepam is a benzodiazepine that may be used for the treatment of panic and seizure disorders.
Product concentration	Lot-specific
Miscellaneous	This product was previously available with a Catalog codes P83-99-21A and LA270. It has been produced using the same conjugation protocol.
Matched pair	100910, 100912, P03-99-21M-P, P03-99-22M-P, HM313 or HM314 (detection)

Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored and used according to the instructions. MedixMAB and MedixAntigens are trademarks of Medix Biochemica and may not be used or reproduced without Medix Biochemica's written permission.

