

Product specification ANTIBODY

2024-07-17

Anti-h NGAL 4203 SPTN-5

Product overview

Catalog number 100580

Specificity Antibody recognizes human neutrophil gelatinase-associated lipocalin

(NGAL)

Description Monoclonal mouse antibody, cultured *in vitro* under conditions free from

animal-derived components.

Product buffer solution 50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN₃ as a preservative

Shelf life and storage 36 months from manufacturing at 2–8 °C

Subclass IgG₁

Analyte description Neutrophil gelatinase-associated lipocalin (NGAL) is a small glycoprotein

expressed in the bone marrow, epithelial tissues associated with

antimicrobial defense, and in the distal tubules and collecting duct of the kidney. Elevated levels of blood and urine NGAL are associated with acute

or chronic renal failure.

Parameters tested on each lot

Product appearance Liquid, may turn slightly opaque during storage

Product concentration 5.0 mg/ml (+/- 10%)

Immunoreactivity 80–120% compared to the reference sample in an FIA test

IEF Profile 6.3–7.3

Purity $\geq 95 \%$

Kinetic parameters

Association rate constant 2.2 x 10⁵ 1/Ms

Dissociation rate constant Does not dissociate under conditions used

Affinity constant Not Applicable (N/A)

Determination method SPR analysis (ProteOn XPR36)

Determination antigen NGAL, Medix Biochemica (Cat 610012)





Product specification ANTIBODY

2024-07-17

Cross-reactivities Not Determined (N/D)

Epitope N/D

Pair recommendations

		DETECTION			
		4202	4203	4204	4205
CAPTURE	4202	-	+	+	+
	4203	+	-	+	+
	4204	+	+	-	+
	4205	+	+	+	-

Following pairs are especially recommended for the below mentioned assays:

IT: 4202 – 4203, 4203 – 4204 and 4203 – 4205

Please note that pair recommendations are based on results obtained by our laboratory. Equally good results may be obtained using other pairs and therefore these recommendations are only indicative.

Platforms tested FIA, IT

Antigens tested Recombinant NGAL antigen, Medix Biochemica 610012 and native NGAL

antigen, Medix Biochemica 342-49.

Product stability TEMPERATURE, TIME RESULT

-70 °C, 21 days OK
-20 °C, 21 days OK
+4 °C, 21 days OK
+35 °C, 21 days OK
+45 °C, 7 days OK

Stability testing is performed in the product buffer to see whether different temperatures affect the antigen binding, charge or composition of the antibody. Please note that the shelf life given on the first page is based on real time stability testing at 2–8 °C in the product buffer.

Miscellaneous -

References -

