

10850 METRO CT. MARYLAND HEIGHTS, MISSOURI, USA 63043 T: 314.968.1091 F: 314.968.9851 WWW.LEEBIO.COM

## SPECIFICATION

# **Neuron Specific Enolase**

**Human Brain** 

#### **GENERAL INFORMATION**

Catalogue Number 430-11 Storage Temperature -20°C

CAS Number 9014-08-8 Storage Note Avoid freeze thaws

E.C. Number 4.2.1.11 Shipping Note Dry Ice

Synonyms NSE

## PRODUCT INFORMATION

Source Human Brain

Form 0.2 micron filtered, frozen solution in 0.01 M sodium phosphate, 0.15 M sodium

chloride and 0.005 M magnesium sulfate, pH 7.4

Appearance Clear, colorless frozen liquid Protein  $\geq 0.5 \text{ mg/mL} (E^{0.1\%}_{280nm} = 1.22)$ 

Purity > 96% (SDS-PAGE)

Purity Note Cluster of bands at 45,000 MW Recertification 5 years from date of manufacture

## HANDLING INFORMATION

Handling Handle as a potentially hazardous substance.

Recovery If small volume, microfuge vial before drawing from container.

Non-infectious disease statement: Negative or non-reactive at the donor level for anti-HIV 1 and 2, anti HCV, and HBsAg by FDA approved methods. However, because no test method can offer complete assurance that infectious agents are absent, this material should be handled at the Biosafety Level 2 (BSL 2) as recommended for any potentially infectious human serum or blood specimen in the CDC/NIH manual "Biosafety in Microbiological and Biomedical Laboratories", 2009. This product is intended for research, laboratory use, or further manufacturing. It is not for administration to humans. We certify this to be true to the best of our knowledge. The statements herein are offered for informational purposes only and are intended to be used solely for your consideration, investigation and verification. The suitability of this product for each application must be determined by each end user.