

Product Data: Influenza B Virus Antibody

Catalog #	HM424
Lot #	B1731
Format	Protein A Purified Monoclonal Antibody
Host/Source	Mouse
Purity	>95% by SDS-PAGE
Immunogen/Antigen	Inactivated Influenza B Virus
Specificity	Influenza B Virus Nucleoprotein
Class/Isotype	IgG2a
Crossreactivity	No reactivity with Influenza A, Adenovirus, and RSV
Manufacture date	3/15/2021
Expiration date	3/14/2023
Concentration	5.1 mg/ml
Supplied In	10 mM Phosphate Buffered Saline, pH 7.4
Preservative	0.05% Sodium Azide
Applications	Lateral Flow
Matched Pair	HM424 (self-pairs)
Storage	2-8°C

Product Data: **Influenza B Virus Antibody**

Note	This product is intended for research and manufacturing uses only. It is not a diagnostic device. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.
-------------	--