

Product Data: Influenza A Virus Antibody

Catalog #	HM418
Lot #	Z2611
Format	Protein A Purified Monoclonal Antibody
Host/Source	Mouse
Purity	>95% by SDS-PAGE
Immunogen/Antigen	Inactivated Influenza A Virus
Specificity	Influenza A Virus Nucleoprotein
Class/Isotope	IgG2a
Crossreactivity	No reactivity with Influenza B, Adenovirus, and RSV
Concentration	4.9 mg/ml
Supplied In	10 mM Phosphate Buffered Saline, pH 7.4
Preservative	0.05% Sodium Azide
Applications	Lateral Flow
Matched Pair	HM419 or HM423 (detection)
Storage	2-8°C
Note	This product is intended for research and manufacturing uses only. It is not a diagnostic device. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.