



211 Wells Street | P.O. Box 489 | North Berwick, ME 03906

P 207.676.7639 | F 207.676.7658 | E info@eastcoastbio.com

www.eastcoastbio.com

Product Data:

Adenovirus Antibody

Catalog #	HM1178
Lot#	B0449
Format	Protein G Purified Monoclonal Antibody
Host/Source	Mouse
Specificity	Hexon Protein
Class/Isotype	IgG2a
Manufacture date	4/30/2020
Expiration date	04/30/2026
Concentration	1.0 mg/ml
Supplied In	25 mM Phosphate Buffered Saline, pH 7.0-7.5
Applications	ELISA, Lateral Flow, IFA, WB
Matched Pair	FG300 (capture)
Storage	2-8°C
Note	This product is intended for research and manufacturing uses only. It is not a diagnostic device. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.