

Product Data: Respiratory Syncytial Virus (RSV) Antibody

Catalog #	HM437
Lot #	Y0807
Format	Protein A Purified Monoclonal Antibody
Host/Source	Mouse
Purity	>95% by SDS-PAGE
Immunogen/Antigen	Purified native RSV
Specificity	Fusion protein of RSV strains A & B
Class/Isotype	IgG2a
Crossreactivity	No reactivity with Influenza A, Influenza B, and Adenoviruses
Manufacture date	3/21/2018
Expiration date	3/20/2022
Concentration	5.0 mg/ml
Supplied In	10 mM Phosphate Buffered Saline, pH 7.4
Preservative	0.05% Sodium Azide
Applications	Lateral Flow
Matched Pair	HM438 or HM439 (detection)
Storage	2-8°C

Product Data:

Respiratory Syncytial Virus (RSV) Antibody

Note	This product is intended for research and manufacturing uses only. It is not a diagnostic device. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.
-------------	--