

EastCoast Bio 211 Wells Street | P.O. Box 489 | North Berwick, ME 03906 P 207.676.7639 | F 207.676.7658 | E info@eastcoastbio.com www.eastcoastbio.com

Product Data: Influenza B Virus Antibody

Catalog #	HM427
Lot #	B2685
Format	Protein A Purified Monoclonal Antibody
Host/Source	Mouse
Purity	>95% by SDS-PAGE
Immunogen/Antigen	Inactivated Influenza B Virus
Specificity	Influenza B Virus Nucleoprotein
Class/Isotype	lgG2b
Crossreactivity	No reactivity with Influenza A, Adenovirus, and RSV
Manufacture date	10/1/2021
Expiration date	10/1/2023
Concentration	5.4 mg/ml
Supplied In	10 mM Phosphate Buffered Saline, pH 7.4
Preservative	0.05% Sodium Azide
Applications	Lateral Flow
Matched Pair	HM428 (detection)
Storage	2-8°C



Product Data: Influenza B Virus Antibody

Note	This product is intended for research and manufacturing uses only. It is not
	a diagnostic device. The user assumes all responsibility for care, custody
	and control of the material, including its disposal, in accordance with all
	regulations.