

Part of **Medix Biochemica Group**

PO Box 489
 North Berwick, ME 03906
 U.S.A.
www.eastcoastbio.com

Product Data: Influenza A Virus Antibody

Product Type	Monoclonal
Catalog #	HM423
Qty/Pack Size	1 mg
Format	Protein A Purified Monoclonal Antibody
Host/Source	Mouse
Purity	>95% by SDS-PAGE
Immunogen/Antigen	Inactivated Influenza A Virus
Specificity	Influenza A Virus Nucleoprotein
Class/Isotype	IgG1
Cross Reactivity	No reactivity with Influenza B, Adenovirus, and RSV
Supplied In	10 mM Phosphate Buffered Saline, pH 7.4
Preservative	0.05% Sodium Azide
Applications	Lateral Flow
Matched Pair	HM418 or HM421 (capture)
Storage	2-8°C
Product Category	Infectious Diseases
Note	This product is intended for research and manufacturing uses only. It is not a diagnostic device. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.

© 2020 EastCoast Bio. All Rights Reserved.