

Part of **Medix Biochemica Group**

PO Box 489
 North Berwick, ME 03906
 U.S.A.
www.eastcoastbio.com

Product Data: Influenza A Virus Antigen, Native

Product Type	Antigen
Catalog #	LA157
Qty/Pack Size	1 ml
Format	Liquid
Host/Source	A/New Caledonia/20/99 (H1N1)
Supplied In	Allantoic Fluid
Preservative	None
Inactivation	Heat treated
Applications	ELISA
Storage	≤-70°C
Product Category	Infectious Diseases
Comments	Although material has been inactivated, as with all biological reagents, it should be treated as potentially infectious.
Note	This product is intended for research and manufacturing uses only. It is not a diagnostic device. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.

© 2020 EastCoast Bio. All Rights Reserved.