

211 Wells Street | P.O. Box 489 | North Berwick, ME 03906

P 207.676.7639 | F 207.676.7658 | E info@eastcoastbio.com

www.eastcoastbio.com

Product Data:

Respiratory Syncytial Virus (RSV) Antibody

Catalag #	LIM424
Catalog #	HM434
Lot #	Y1294
Format	Protein A Purified Monoclonal Antibody
Host/Source	Mouse
Purity	>95% by SDS-PAGE
Immunogen/Antigen	Purified native RSV
Specificity	Fusion protein of RSV strains A & B
Class/Isotope	IgG1
Crossreactivity	No reactivity with Influenza A, Influenza B, and Adenoviruses
Concentration	5.0 mg/ml
Supplied In	10 mM Phosphate Buffered Saline, pH 7.4
Preservative	0.05% Sodium Azide
Applications	Lateral Flow
Matched Pair	HM433 (capture)
Storage	2-8°C
Note	This product is intended for research and manufacturing uses only. It is not a diagnostic device. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.