

## Product Data: **Adenovirus Antibody**

<b>Catalog #</b>	HM056
<b>Lot #</b>	A0751
<b>Format</b>	Purified Monoclonal Antibody
<b>Host/Source</b>	Mouse
<b>Purity</b>	>95%
<b>Immunogen/Antigen</b>	Purified human adenovirus
<b>Specificity</b>	Adenovirus
<b>Class/Isotype</b>	IgG2a
<b>Manufacture date</b>	12/1/2017
<b>Expiration date</b>	12/1/2022
<b>Concentration</b>	6.91 mg/ml
<b>Supplied In</b>	10 mM Phosphate Buffered Saline, pH 7.0
<b>Preservative</b>	0.02% Sodium Azide
<b>Applications</b>	Lateral Flow
<b>Matched Pair</b>	HM057 (detection)
<b>Storage</b>	Short Term ( $\leq$ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

## Product Data:

# Adenovirus Antibody

<b>Note</b>	This product is intended for research and manufacturing uses only. It is not a diagnostic device. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.
-------------	--