

Product Data:

**Respiratory Syncytial Virus (RSV-fusion) Antigen, Recombinant  
>90%**

<b>Catalog #</b>	LA512
<b>Lot #</b>	C0524
<b>Format</b>	Purified Recombinant Antigen (His-tag)
<b>Host/Source</b>	Mammalian Cell Expression System
<b>Purity</b>	>90%
<b>Manufacture date</b>	2/9/2022
<b>Expiration date</b>	08/31/2022
<b>Concentration</b>	1.71 mg/ml
<b>Supplied In</b>	Acetate Buffer with salts and other additives, pH 5.0
<b>Preservative</b>	0.09% Sodium Azide
<b>Applications</b>	ELISA, Lateral Flow
<b>Matched Pair</b>	HM947
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Note</b>	This product is intended for research and manufacturing uses only. It is not a diagnostic device. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.