

## Product Data:

# Respiratory Syncytial Virus (RSV) Antigen, Native

<b>Catalog #</b>	FC517
<b>Lot #</b>	X1021
<b>Format</b>	Liquid
<b>Host/Source</b>	Hep-2 cells/Strain Long
<b>Manufacture date</b>	3/2/2017
<b>Expiration date</b>	3/2/2022
<b>Concentration</b>	1.0 mg/ml
<b>Supplied In</b>	50 mM Glycine Buffer with 0.05% CHAPS and 0.004% SDS
<b>Preservative</b>	None
<b>Inactivation</b>	Photochemical and Detergent
<b>Applications</b>	ELISA
<b>Storage</b>	≤-70°C
<b>Comments</b>	Although inactivated, material should be treated as potentially infectious.
<b>Note</b>	This product is intended for research and manufacturing uses only. It is not a diagnostic device. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.