

## Product Data: Hemoglobin (Hb) Antibody

<b>Catalog #</b>	HM258
<b>Lot #</b>	Z2813
<b>Format</b>	Ion Exchange Purified Monoclonal Antibody
<b>Host/Source</b>	Mouse
<b>Purity</b>	>95%
<b>Specificity</b>	Human hemoglobin
<b>Class/Isotope</b>	IgG1
<b>Crossreactivity</b>	No reactivity with hemoglobin from pig, cow, goat, sheep, horse, chicken, rabbit
<b>Manufacture date</b>	08/01/2019
<b>Expiration date</b>	2/1/2022
<b>Concentration</b>	22.68 mg/ml
<b>Supplied In</b>	10 mM Phosphate Buffered Saline, pH 8.0 with 50 mM NaCl
<b>Preservative</b>	0.1% Sodium Azide
<b>Applications</b>	ELISA, Lateral Flow
<b>Matched Pair</b>	HM257 (detection)
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Comments</b>	Sensitivity of ≤ 30 ng/ml

## Product Data: **Hemoglobin (Hb) Antibody**

<b>Note</b>	This product is intended for research and manufacturing uses only. It is not a diagnostic device. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.
-------------	--